NEUROGUIDE SIGNAL GENERATOR MANUAL AND TUTORIAL

Copyright © 2003-2006 Applied Neuroscience, Inc. (Sine Wave segments were selected for illustrative purposes only)

Introduction:

The signal generator is used to calibrate and test the digital signal processing properties of NeuroGuide as well as to serve as an educational program by which the principles of digital signal analyses can be learned and explored. Concepts such as frequency, time, phase delays, noise, amplitude and coherence can be tested and evaluated.

EEG data can be simulated by approximating the selected mixtures of signals to match the signal parameters and scalp locations of the EEG.

TABLE OF CONTENTS

- Step #1 Launch NeuroGuide and click File>Open>Signal Generation
- Step #2 <u>Use Mouse to Select EEG Channels</u>, Sine Wave Frequencies and Amplitudes (uV) and Phase Delays (degrees) and "Noise" (% S/N ratio)
- Step #3 Simulate EEG 'Spindles' using the Pulse generate option
- Step # 4- Click OK, then Click Edit>Select All to view FFT results
- Step # 5 <u>Click File>Save As</u> to save the signals in NeuroGuide or Lexicor format (*.ng or *.dat).
- Step #6 <u>Example Tutorial</u> of Replicating Peer Reviewed Publication: Gomez and Thatcher "Frequency domain equivalence between potentials and currents using LORETA." <u>Int. J. of Neuroscience</u>, 107: 161-171, 2001.
- <u>Appendix A</u> LORETA
- **<u>Appendix B</u>** Mathematics of Gomez and Thatcher, 2001
- <u>Appendix C</u> References

Step #1: Launch NeuroGuide and click File>Open>Signal Generation

Return to Top

Notar of Palitie		(6).
the Edk Inter Blan	lage Beport Analysis Statistic	a Wedee Seb
Qperi •	NeuroGuide Obten	
	HESA +	
Lago ine	Bolens	
	Brantfaller *	
	Davel.	
	Çedicen 🕐	
	Devned	
Est .	Leskie ·	
	ignizer	
100	Among Among and	
	Americana.	
NEW AND AND A STREET AND A ST	ignorificat	
337	Subips	
	boCorp Export File	
· · · · · · · · · · · · · · · · · · ·	GP	
ra calle	Section.	
	Sector Sector	
	Contractor and the second	
	Source product	
about the		
Trooling		
5		

Step #2 - Use Mouse to Select EEG Channels, Waveforms (sine, sawthooth, square and pulses), Frequencies (1 to 30 Hz), Amplitudes (uV), Phase Delays (degrees) and "Noise" (% S/N ratio)



2a - Click a channel to select a location in Lexicor format (e.g., O1), then double click a Frequency (e.g., 10 Hz), then double click Amplitude (uV) and type in a value (e.g., 100 uV).

ignal	Signal Pa	rameters by Frequen	cy.	Signal Waveform
FP1 FP2		Amplitude (uV)	Phase Shift (Deg)	
F3	0 Hz	0.00	0.00	
F4	1 Hz	0.00	0.00	
ČĂ	2 Hz	0.00	0.00	
P3 P4	3 Hz	0.00	0.00	
DÎ	4 Hz	0.00	0.00	
02 57	5Hz	0.00	0.00	1114114111411141114
8	6 Hz	0.00	0.00	
13	7 Hz	0.00	0.00	1 1 1 1 1 1 1 1 1 1 1 1 1
15	8 Hz	0.00	0.00	
16	9 Hz	0.00	0.00	********
-2 Cz	10 Hz	100.00	0.00	11
Pz	11 Hz	0.00	0.00	
	12Hz	0.00	0.00	
ignal Noise (uV)	13 Hz	0.00	0.00	QK Cancel
0.00	14 Hz	0.00	0.00	
90.00	15 Hz	0.00	0.00	

2b – Mix sine waves in by double clicking the amplitude of a different frequency, e.g., 5 Hz and type 50 uV.

iignal	Signal Pa	rameters by Frequen	2	Signal Waveform
FP1 FP2		Amplitude (uV)	Phase Shift (Deg)	
F3	0 Hz	0.00	0.00	
-4	1 Hz	0.00	0.00	1.1.1.1.1.1.
24	2 Hz	0.00	0.00	
3	3 Hz	0.00	0.00	11111111111111111111
n I	4 Hz	0.00	0.00	
12	5 Hz	50.00	0.00	
8	6 Hz	0.00	0.00	V (
3	7 Hz	0.00	0.00	
5	8 Hz	0.00	0.00	
6	9 Hz	0.00	0.00	V V V
2	10 Hz	100.00	0.00	
2	11 Hz	0.00	0.00	
	12 Hz	0.00	0.00	· · · · · · · · · · · · · · · · · · ·
gnal Noise (uV)	13 Hz	0.00	0.00	QK Cancel
0.00	14 Hz	0.00	0.00	
0.00	15 Hz	0.00	0.00	

2c – Shift the Phase of the 5 Hz signal by double clicking "Phase Shift (Deg)" at 5 Hz and type 30.

ignal	Signal Pa	rameters by Frequen	cy.	Signal Waveform
FP1 FP2		Amplitude (uV)	Phase Shift (Deg)	
3	0 Hz	0.00	0.00	A A A A A
4	1 Hz	0.00	0.00	B A B A B A B A B A
Ä	2 Hz	0.00	0.00	
3	3 Hz	0.00	0.00	
ที่ 1	4 Hz	0.00	0.00	
02 7	5Hz	50.00	30.00	
8	6 Hz	0.00	0.00	
13	7 Hz	0.00	0.00	
5	8 Hz	0.00	0.00	
6	9 Hz	0.00	0.00	V V V V
2	10 Hz	100.00	8.00	
Pz .	11 Hz	0.00	0.00	
	12Hz	0.00	0.00	
ignal Noise (uV)	13 Hz	0.00	0.00	QK Cancel
0.00	14 Hz	0.00	0.00	
Jacob.	15 Hz	0.00	0.00	

2d – Add "Noise" to the 5 Hz signal by double clicking the window below "Signal Noise (uV)" and type 100. This adds 100 microvolts of noise to the 5 Hz signal located at O1.

ignal	Signal Parameters by Frequency			Signal Waveform
FP1 FP2		Amplitude (uV)	Phase Shift (Deg)	I I
F3	0 Hz	0.00	0.00	
F4	1 Hz	0.00	0.00	
C4	2Hz	0.00	0.00	- 附前 方面 方面 防病
P3 P4	3 Hz	0.00	0.00	THE ALCONDATION TO DEPART
01	4 Hz	0.00	0.00	- / 1 x// +
02 57	5 Hz	50.00	30.00	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
F8	6 Hz	0.00	0.00	
13	7 Hz	0.00	0.00	
15	8 Hz	0.00	0.00	1 1 N 1
TE	9 Hz	0.00	0.00	
rz Cz	10 Hz	100.00	0.00	172
Pz	11 Hz	0.00	0.00	
	12Hz	0.00	0.00	0
ignal Noise (uV)	13 Hz	0.00	0.00	QK Cancel
-	14 Hz	0.00	0.00	
(Jacobia	15 Hz	0.00	0.00	

2e - Repeat Steps 2a to 2d for each channel with or without adding phase delays and/or noise or multiple frequencies, etc. Unselect any value by double clicking in the appropriate box and set the value = 0. The Channel is the primary selection and then the amplitude, frequency or mixtures of frequencies and phases and noise are the secondary selections.

Step #3 - Simulate EEG "Spindles" by selecting pulse and then select the frequency and amplitude of the intra-pulse structure, duration and inter-pulse intervals of the simulated spindles.

3a - Generate Pulses of different durations and inter-pulse intervals. Use this tool to simulate EEG "Spindles".





3b - Simulate any EEG by comparing the auto and cross-spectral values and then entering these values into the appropriate channels and appropriate parameter selection locations. Use the Signal Generator feature of NeuroGuide to learn about digital signal processing in general as well as various analytical programs including LORETA and other inverse solutions.

Step # 4- Click OK, then Click Edit>Select All to view FFT results



Step #5 - Click File>Save As to save the signals in NeuroGuide or Lexicor format (*.ng or *.dat).

Return to Top

Step – 5a - Follow the NeuroGuide Manual Instructions (step #6) to save as NeuroGuide (*.ng) or Lexicor (*.dat) files.

Step – 5b - Follow the NeuroGuide Manual Instructions (step # 6) to save Power Spectra and Cross-Spectra (Step # 6) and to Export to LORETA (Step # 11 in the NeuroGuide Manual).

Step #6 - Example Tutorial by replicating the publication: Gomez and Thatcher "Frequency domain equivalence between potentials and currents using LORETA." <u>Int. J. of</u> <u>Neuroscience</u>, 107: 161-171, 2001.

Signal	Signal Pa	rameters by Frequen	2	Signal Waveform
FP1 FP2		Amplitude (uV)	Phase Shift (Deg)	
F3 F4	7 Hz	0.00	0.00	
ci	8 Hz	0.00	0.00	
C4	9 Hz	0.00	0.00	
P4	10 Hz	50.00	0.00	
01	11 Hz	0.00	0.00	
F7	12 Hz	0.00	0.00	
F8	13 Hz	0.00	0.00	
T4	14 Hz	0.00	0.00	
15	15 Hz	0.00	0.00	
16 F2	16 Hz	0.00	0.00	
Cz	17 Hz	0.00	0.00	
P2	18 Hz	0.00	0.00	
	19 Hz	0.00	0.00	1
Signal Noise (uV)	20 Hz	25.00	0.00	QK Cancel
0.00	21 Hz	0.00	0.00	
Sec.	22 Hz	0.00	0.00	

5a- Select O1 at 10 Hz at 50 uV and 20 Hz at 25 uV

$5b-Select\ O2$ at 8 Hz 50 uV and 16 Hz at 25 uV



$5c-Select\ F7$ at 10 Hz 25 uV and 20 Hz at 50 uV

Signal	Signal Pa	rameters by Frequen	cy.	Signal Waveform
FP1 FP2		Amplitude (uV)	Phase Shift (Deg)	LAND LAND
F3 F4	7 Hz	0.00	0.00	
	8 Hz	0.00	0.00	1.1.1.1. D. M. (. 1. 1. 1. 1. 1.
C4	9 Hz	0.00	0.00	一封作任意的准认作只是当作到多的准法推出。
P4	10 Hz	25.00	0.00	- 114 114 114 114 114 114 114 114 114 11
01	11 Hz	0.00	0.00	
6/2 F7	12 Hz	0.00	0.00	
F8	13 Hz	0.00	0.00	
T4	14 Hz	0.00	0.00	
T5	15 Hz	0.00	0.00	
16 F2	16 Hz	0.00	0.00	
Cz	17 Hz	0.00	0.00	
P2	18 Hz	0.00	0.00	
	19 Hz	0.00	0.00	1
Signal Noise (uV)	20 Hz	50.00	0.00	QK Çancel
0.00	21 Hz	0.00	0.00	
850	22 Hz	0.00	0.00	

5d – Select F8 at 8 Hz 25 uV and 16 Hz at 16 uV

Signal	Signal Parameters by Frequency			Signal Waveform
FP1 FP2		Amplitude (uV)	Phase Shift (Deg)	
F3 F4	7 Hz	0.00	0.00	
iii ii	8 Hz	25.00	0.00	1. 1. 1. 1. 1. 1. 1. 1. 1. 1. 1. 1. 1. 1
C4 P3	9 Hz	0.00	0.00	- 118 118 138 118 118 118 118 118
P4	10 Hz	0.00	0.00	- 作我们我们我们我们我们我们我们
01	11 Hz	0.00	0.00	
F7	12 Hz	0.00	0.00	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
F8	13 Hz	0.00	0.00	A [] A [] A [] A [] A [] A [] A [
T4	14 Hz	0.00	0.00	
15	15 Hz	0.00	0.00	
Fz	16 Hz	50.00	0.00	
Cz	17 Hz	0.00	0.00	
P2	18 Hz	0.00	0.00	
	19 Hz	0.00	0.00	
Signal Noise (uV)	20 Hz	0.00	0.00	QK Çancel
0.00	21 Hz	0.00	0.00	
	22 Hz	0.00	0.00	

5e – Click OK and then click File>Save As NeuroGuide (*.ng) or Lexicor (*.dat) for purposes of further analysis. For example, save the power spectra and/or export the digital time series to LORETA.



APPENDIX – A

Return to Top

Gomez and Thatcher, 2001 used the Key Institute mathematical equations to validate LORETA and cross-validated their math by comparing the forward solution and the inverse solution using MRI 3-D voxel locations and the surface scalp EEG currents and potentials (Based on the Reciprocity Theorum, Helmholtz, 1853). The results of the Gomez and Thatcher, 2001 study is also consistent with Tesche, C. and Kajola, M. "A comparison of the localization of spontaneous neuromagnetic activity in the frequency and time domains." <u>Electroencephalography and Clin. Neurophysiology</u>, 87(6): 408-416, 1993.

One can test the facts and the science of LORETA for themselves using the NeuroGuide signal generator and the Gomez and Thatcher, 2001 frequencies and locations which is only one of several tools available to test LORETA (see Appendix B and C) not to mention the mathematical concepts of linearity between frequency and time and between electrical potentials and currents (Helmholtz, 1853 physics of "Reciprocity" and the "Lead Field", Malmivuo and Plonsey, 1995).

It makes no difference whether one exports signals in the time domain or in the frequency domain (as demonstrated in the Gomez & Thatcher, 2001 and the Tesche et al, 1993 publications as well as by mathematical simulation in step # 5). Caution must be exercised when using LORETA to be sure to physiologically validate using the surface linked ears, average reference and current source density data. This is not to indicate that LORETA is not a valid mathematical and scientific methodology, to the contrary, it is an important contribution. We are emphasizing the fact that LORETA is valid when used by competent professionals who take the time to validate the source solutions by evaluating the surface EEG distributions and physiology in order to guard against localization error. For example, scalp recorded EEG with large amplitude alpha in O1 and O2 should appear as high current density in areas 17, 18 & 19 in LORETA.

APPENDIX – B Mathematics and Results of Gomez & Thatcher, 2001

Return to Top

Note: There are three instances when multiplication of matrices is communitative: 1- by a null matrix, 2by an identity matrix and, 3- multiplication by a scalar. The equation below is a valid equality when using a scalar as we do.

$$\lambda A = \{\lambda a_{ij}\} = \{a_{ij}\lambda\} = A\lambda$$
 Eq. 1

We apply this community property in the following manner. For S = KJ, where K is the lead field matrix, J = current and S = the sensitivity of the sensors (dependingon the model used and the conductivity, etc.). S is an N x W matrix for the scalp potentials (EEG/MEG),where N is the number of sensors and W is the number of time samples. J is a 3M x W matrix, where Mis the number of sources and W is the same time samples as for S. Then the inverse solution is a linearcombination of the signal S in the sensors

$$\mathbf{J} = \mathbf{T} \cdot \mathbf{S} \qquad \qquad \mathbf{Eq. 2}$$

Where T is some generalized inverse of K, where the minimum norm solution is

 K' is the transpose of K, and · represents matrix multiplication and pinv(X) is the Moore-Penrose pseudoinverse (Menke, 1984). Pascual-Marqui et al, 1994 use the mathematical method that they refer to as "Low-Resolution Computed Tomography" (LORETA) to add physiological foundations and to avoid the minimum norms's problems in localizing deep sources by using the Laplacian Operator B and W as a weighting matrix. The LORETA equation is

 $T = \{pinv(WB'BW)\}K'[pinv(K inv(WB'BW)K']$ Eq. 4

The critical factor in these considerations is that the real number FFT computed by the cross-spectrum (Hermitian matrix as a scalar real number) as represented in equation 1 is a linear operator such that for any inverse solution of the form in equation 3 is equivalent to:

$$FFT(J) = FFT [T \cdot S] = T \cdot FFT[S]$$
Eq. 5

Equation 5 is the formula that Gomez and Thatcher (2001) used. Gomez and Thatcher (2001) simulated the linear equivalence by a combination of sine waves and confirmed the linearity of equation 5 as any one can do by using the NeuroGuide signal generator as described in step # 5 for oneself.

Figure 1 – From Gomez & Thatcher, 2001. This is the three-shell spherical model of the head used to simulate LORETA. Four electrodes (F7, F8, O1, O2) and the reference electrode (A1) are indicated by black rectangles. The coordinates of the electrodes are according to the best-fitting sphere relative to cortical anatomy (Towel et al., 1993). The peaks of beta (for F7 and F8) and alpha activity (for O1 and O2) are indicated in parenthesis. Eight sources (1 to 3) indicated by black circles were located in the interior of the sphere to represent the equivalent current sources such as in the gray matter.



Figure 2 – From Gomez & Thatcher, 2001. Power spectrum of the signals used to simulate LORETA. The spectrum of the signals in the scalp electrodes is shown on the left (amplitude of beta is higher in the anterior regions, alpha amplitude is higher in the posterior regions and a frequency shift toward the lower frequencies in the right hemisphere). The center and right columns are the spectra of the current sources nearest to the electrodes J1, J3, J5 and J7 after calculating the inverse solution. Each source has three components x, y and z. The y-axis of the electrodes is uv^2 /cycle/sec for the electrodes and uA/cm^2) ²/cycle/sec for current density at each source location. The x-axis is frequency in Hz in all cases. This simulation confirms the mathematical statements and demonstrates a frequency domain equivalence between the spectra of electrical potentials at the scalp and the spectra of currents in the interior of the head model.



APPENDIX – C - REFERENCES

Return to Top

Baillet, S., et al, "Evaluation of inverse methods and head models for EEG source localization using a human skull phantom". Physics in Medicine and Biology, 46: 77-96, 2001.

Baillet, S. and Garnero, L. "A Bayesian approach to introducing anatomo-functional priors in the EEG, MEG inverse problem". IEEE Trans. Biomed. Eng. 44: 374-375, 1997.

Casper, et al. "Evaluation of inverse methods and head models for EEG source localization using a

human skull phantom" at: http://sipi.usc.edu/~silvin/docs/inversecasperthese2.pdf

Gomez, J. F. and Thatcher, R.W. "Frequency domain equivalence between potentials and currents using LORETA." Int. J. of Neuroscience, 107: 161-171, 2001.

Helmholtz, HLF, Ann. Physik and Chemie 89: 211-233, 354-377, 1853 (see also "Helmholtz's Treatise on Physiological Optics" by Hermann Von Helmholtz, edited by J. P. Southal, Thoemmes, Press, 2000, ISBN 1855068311).

Malmivuo, J. and Plonsey, R. "Bioelectromagnetism", Oxford University Press, 1995.

Menke, W. "Geophysical Data Analysis: Discrete Inverse Theory." Orlando: Academ, ic Press, 1984.

Hämäläinen, *M.* "Discrete and distributed source estimates",*ISBET Newsletter* Edited by W. Skrandies, Giessen, Germany, No 6 / December 1995.

Pascual-Marqui, R. D. Review of methods for solving the EEG inverse problem. Inter. J. of Bioelectromagnetism, 1:75-86, 1999.

Tesche, C. and Kajola, M. "A comparison of the localization of spontaneous neuromagnetic activity in the frequency and time domains." <u>Electroencephalography and Clin. Neurophysiology</u>, 87(6): 408-416, 1993.

Towel, V. et al., The spatial location of EEG electrodes: locating the best-fitting sphere relative to cortical anatomy. EEG & Clin. Nerurophysiol., 103: 9 – 15, 1993.